

## **Micro-PET Hot-Rod Phantom**

## Measure Spatial resolution of small animal micro-PET and SPECT systems.

The QRM-MicroPET-HotRod Phantom is a perfect tool to evaluate spatial resolution of PET and SPECT systems designed for small animal imaging.

The phantom is designed for a FOV smaller then 40 mm and is made completely of PMMA. The phantom consists of a fillable cylinder for the isotope.

Two spacers can be used to align the resolution insert in the center of the phantom. The spacers can be exchanged by two full body PMMA inserts in order to avoid an excessively high signal of the isotope inside the phantom

The top of the phantom houses a bubble reservoir to avoid air bubbles inside the scan area if placed in horizontal position in the scanner.

## **Specifications**

phantom
---------

Base material	PMMA
Diameter	35 mm
Height	70 mm

## HotRod insert

Diameter 29 mn	n
Height 12 mm	n
Hole pattern 0.6, 0.8, 1.0, 1.2, 1.5, 2.0 mm	n
Center to center distance 2 x diamet	e

Ø	Ø	qy	ctc dist.
0.6 mm	0.024 in.	43	1.2 mm
0.8 mm	0.031 in.	26	1.6 mm
1.0 mm	0.039 in.	19	2.0 mm
1.2 mm	0.047 in.	13	2.4 mm
1.5 mm	0.059 in.	8	3.0 mm
2.0 mm	0.079 in.	6	4.0 mm





MicroPET-HotRod Phantom



parts included

