Multi-Purpose Phantom Model ATS 539



The Model 539 Multipurpose phantom is an easy, comprehensive means of evaluating imaging systems over the full range of clinical imaging frequencies (2 MHz to 18 MHz).

The phantom is designed with a combination of monofilament line targets for distance measurements and tissue mimicking target structures of varying sizes and contrasts. Cystic-like target structures are positioned in-line vertically, thereby permitting an entire target group to be displayed in one view. Due to the acoustic similarity of the background material and the target structures, artifacts caused by distortion, shadowing or enhancement have been eliminated. Six gray scale targets ranging in contrast from +15 to -15 dB are provided to evaluate the system's displayed dynamic range and gray scale processing performance. All ATS urethane phantoms are guaranteed for the useful life of the phantom, defined as 10 years.


Key Tests with Model ATS 539

- Uniformity
- Depth of Penetration
- Beam Profile/ Focal Zone/ Lateral Response
 Width
- Vertical Distance Measurement
- Horizontal Distance Measurement
- Axial and Lateral Resolution
- Elevational Resolution
- Contrast Resolution
- Grayscale Contrast Sensitivity
- Dead Zone Assessment



深圳市展业达鸿科技有限公司

吕先生:15920060912(微信同号) 0755-22934005(座机) 地址:深圳市福田区八卦二路八卦岭工业区615 栋419 邮箱:hongqi@thingstet.com 网址:www.thingstest.com



SPECIFICATIONS

DIMENSIONS	23.4 x 20.5 x 9.5cm (9" x 8" x 4")
PHANTOM WEIGHT	7 lbs. (3.1 kg)
HOUSING MATERIALS	PVC
WALL THICKNESS	1.0 cm
SCAN SURFACES	4
SCAN SURFACE DIMENSIONS	17.5x7.5 cm 14.0x7.5 cm
TISSUE-MIMICKING MATERIAL	Urethane Rubber

URETHANE PROPERTIES

< -40°C Freezing point: Melting point: Above 100° C Speed of Sound: 1450 m/s at 23° Attenuation Coefficient: 0.5 dB/cm/MHz (measured at 3.5 MHz) Above 100° C

LINE TARGETS

Material: Diamter:

VERTICAL DISTANCE GROUP

Number of targets: Depth range: 17 Spacing: 1.0 cm

HORIZONTAL DISTANCE GROUP

Number of groups: Depths: 2 Scanning Surface #1: 5.0 cm Scanning Surface #3: 5.0 cm Number of Targets: 10 2.0 cm Spacing:

DEAD ZONE GROUP

Number of targets: Depth range: Spacing:

Monofilament nylon 0.12 mm

1.0 - 18.0 cm

2.0-10.0mm 1.0 mm

©2013 Computerized Imaging Reference Systems Inc. All rights reserved. Specifications subject to change without notice. Publication: ATS 539 DS 032420

MODEL ATS 539 INCLUDES

QTY	COMPONENT DESCRIPTION
1	Multi-Purpose Phantom
-	10-Year Warranty
-	User Guide
-	Certificate of Compliance

AXIAL-LATERAL RESOLUTION GROUPS

Number of Targets: 11 Lateral Displacement: 1.0 mm Spacing: 4, 3, 2, 1, 0.5 & 0.25 mm Depth: 7.0, 11.0, 4.0, 16.0 cm

ANECHOIC STEPPED CYLINDERS

Type: Non-echogenic, cylindrical Number of Groups: 5

Group 2, 3, 4: Number of targets: 2.0, 3.0, 4.0 mm Diameter of targets: 1.0 - 17.0 cm 1.0 cm

6.0, 8.0 mm

1.0 - 17.0 cm 2.0 cm

Group 6, 8: Number of targets: Diameter of targets: Depth of Targets: Spacing:

Depth of Targets: Spacing:

GRAY SCALE TARGETS

Type: Number of Targets: Spacing: Depth: Contras

Echogenic, Cylindrical 6 2.0 cm (center to center) 4.0 cm

+15+6, +3-3-6, -15 dB with respect to back 司 深圳市展业达鸿科技有限公司

TERF

旨先生中ut\$592000609€12€(織的宮南市,号c)has 0€7550rii22984005(座机) to (ISO) 13485:2016. 业地验深圳市福田区八卦二路八卦岭工业区615栋419

邮箱:hongqi@thingstet.com 网址:www.thingstest.com